# Instruton E

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enzymatic cleaning concen-  
trate for chemothermic and  
manual reprocessing of  
endoscopes and general  
instrument cleaning

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## Range of application:

Instruton E is an enzymatic instrument cleaner for the chemothermic and manual reprocessing of endoscopes and medical instruments.

## Special characteristics:

- very effective and efficient cleaning
- pH neutral
- biodegradable
- good material compatibility
- very efficient use

## Application:

- chemothermic:  
The dosing of 0.5 – 1.0 % is regulated via the pumps inside endoscope cleaning and disinfecting machines. Recommended water temperature: 30 – 40 °C
- manual:  
Choose a dosage between 0.5 – 3 % according to the degree of soiling. Residue free within 2 – 3 minutes. Optimum water temperature: 30 – 40 °C

## Ingredients:

protease, amylase, surfactants, solubiliser, corrosion inhibitor

## Physical and chemical characteristics:

Appearance: clear, slightly yellow liquid  
Scent: characteristic  
pH value: approx. 6 at 18 – 28 °C  
Density: approx. 1.03 g/cm<sup>3</sup> at 20 °C

## Forms of delivery:

- plastic container of 2 L or 5 L

This product is manufactured according to the international Quality Management System EN ISO 9001.

Further information is provided by your expert advisor upon request or on our homepage.

# ANTISEPTICA

**ANTISEPTICA chem.-pharm. Produkte GmbH**  
D-50259 Pulheim, Carl-Friedrich-Gauß-Straße 7  
Tel.: +49-2234-984 66-0, Fax: +49-2234-984 66-11  
E-Mail: info@antiseptica.com, www.antiseptica.com

Quality Management  
Medical Device  
Directive 93/42/EEC

Quality Management  
We are certified  
Voluntary participation in regular  
monitoring according to ISO 9001:2008

Quality Management  
Voluntary participation in regular  
monitoring according to  
EN ISO 13485:2003 + AC:2007

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**ANTISEPTICA Austria GmbH**