

# Aero-flex

The next nebulizer generation

**AERO-FLEX** 

[www.aero-flex.eu](http://www.aero-flex.eu)  
www.aero-flex.eu

# Aero-Flex

The next nebulizer generation

- Small
- Silent
- Portable
- Easy to use
- Light weight



# Small, light, Portable

Always with you

Easy to manipulate

Practical to carry



9 cm



8,5 cm



120 gr

# Silent

A ultrasonic nebulizer

For inhalations

Quiet

Peaceful

Un-disturbing



Z  
Z  
Z  
Z

**AERO-FLEX**

[www.aero-flex.eu](http://www.aero-flex.eu)  
[www.aero-flex.eu](http://www.aero-flex.eu)



# Easy to use

## One button

To start

To stop

To adjust in accordance  
the 3 intensities



# Technical characteristics







- Maximum consumption: 5W
- Ultrasonic frequency: 2,5 MHZ
- Particles size: between 1 and 5 micron
- Three intensities of mist
  - In low intensity: 0,3ml / min
  - In medium intensity: 0,5ml / min
  - In high intensity: 0,75ml / min
- Made in EEC (European countries)
- Micro-controlled
- Memory of the last selected intensity
- Automatically turns-off when the mist stops



**CE 0197**  
*Notified body T.U.V.*

# Our box



- ☑ 1  Aero-flex
- ☑ 1  car adapter
- ☑ 1  AC / DC adapter
- ☑ 1  cup measurer
- ☑ 1  infant mask
- ☑ 1  adult mask
- ☑ 1 instruction for use
- ☑ 1 year warranty

# Our brand

Aero-flex



Aero-flex classic - AE 1000



Aero-flex customized - AE 1001



Aero-flex customized - AE 1002



Aero-flex customized - AE 1003

Accessories



Kit Trachea - AE 1140



AC/DC adapter - AE 1200



car adapter - AE 1210



Cup measurer - AE 1130



Adult mask - AE 1110



infant mask - AE 1120