




Aukey

Ultra-Portable Video Magnifier


*Travel light,
Color your life!*



Key Features

- ▶ Color, B/W and Reverse images
- ▶ 2 levels, 1.5X to 17X magnification range
- ▶ User-friendly tactile buttons
- ▶ 3.5" Wide-screen TFT LCD
- ▶ Freeze frame
- ▶ Only 3 buttons, very easy to use
- ▶ Aumed Anti-Blur™ Technology
- ▶ Rechargeable Li-Ion battery with 2500mAh capacity, approx. 2.5 hours playback time
- ▶ Lightweight 120g (4.2 Oz)
- ▶ Dimensions 105×72×17mm
- ▶ Two-year manufacturer's warranty

Aukey Specifications



Zooming
Magnification range:
1.5x-17x(2 magnification levels)



Image freezing function available



3 color modes
(True color;B&W; W&B)



Battery working time
2.5hours



Screen size
3.5" TFT LCD

Aumed Group Corp.

www.aumed.cn
www.aumedgroup.com



FC CE