

Aumax-LCD

Multi-purpose Desktop Video Magnifier

The only multi-purpose low vision viewer designed with Aumed Anti-Blur™ Technology. Ideal for both left and right-handed students, with dual control panels on the revolving camera. This advanced video magnifier can provide student's desktop and teacher's "chalk board" magnification.

Camera head slides fully from left to right side of LCD Magnifier.



Camera Head



Camera Pannel



Control Pad

Key Features

- ▶ 19" wide screen LCD
- ▶ 3.5X - 79X magnification
- ▶ 8 view display modes
- ▶ Auto focus
- ▶ Seek & find capability
- ▶ Dual VGA outputs
- ▶ Travel case (optional)
- ▶ Remote control pad

Additional functions with the control pad:

- ▶ Freeze frame
- ▶ Focus lock
- ▶ White balance
- ▶ Line marker (interval adjustable)
- ▶ Left/right hand operation
- ▶ 3 in 1 function: Reading, Distance Viewing, Self Viewing



Aumax-LCD Specifications



Seek and Auto-find function



Weight approx.: 14KG
Size: 55*44*20cm (W*H*D)
(19inch TFT LCD Monitor)



Books, photos, catalogues,
newspapers, magazines etc



Chequebooks
Crosswords
Letters



Easy operation with Line-in Control Pad:
*Freeze frame
*Focus lock
*White Balance
*Switch of auto & manual focus
*Line marker & mask function
*Left/right hand operation mode



Optional
* X-Y Reading Table
* Carrying Case (specially designed for
Aumax-LCD with impact resistance,
shock-proof, portable mobility functions)



Zooming
Magnification range
3.5x-79x



Mode
8 viewing mode



Auto focus function



Screen size
19" TFT LCD

Aumed Group Corp.

www.aumed.cn
www.aumedgroup.com



Aumed Authorized Distributor

