

Aumax-VGA

Multi-purpose Desktop Video Magnifier



- ▶ With the Aumed Anti-Blur™ Technology, Aumax-VGA provides clear, smooth and bright images.
- ▶ All-in-one machine for reading, writing, hand working and self-viewing etc.
- ▶ When connected to a digital projector or PC, Aumax-VGA is able to satisfy both personal and business requirements.
- ▶ Optional travel case with solid design, ideal for mobile applications.
- ▶ The Dual control panels on the top & bottom of the camera make it easy to be operated by either the left or right-handed users.



Control Pad



Camera Head



Camera Panel

Key Features

- ▶ 3.5X - 79X magnification
- ▶ 8 viewing display modes
- ▶ Auto focus
- ▶ Seek & find capability
- ▶ Dual VGA outputs
- ▶ Travel case (optional)
- ▶ Remote control pad

Additional functions with the control pad:

- ▶ Freeze frame
- ▶ Focus lock
- ▶ White balance
- ▶ Auto-Manual focus
- ▶ Line marker (interval adjustable)
- ▶ Left/right hand operation

Aumax-VGA

Specifications



Zooming
Magnification range:
3.5x-79x



Mode
8 viewing modes



Auto focus function



Display:
Without Monitor with Aumax-VGA
Connection with display devices:
by VGA(2 VGA output)



Seek and Auto-find function



Weight: Approx. 2.85KG
Size: 38*22*47cm (W*H*D)



Books, photos, catalogues,
newspapers, magazines etc



Chequebooks
Crosswords
Letters etc.



Easy operation with Line-in Control Pad:
*Freeze frame
*Focus lock
*White Balance
*Switch of auto & manual focus
*Line marker & mask function
*Left/right hand operation mode



Optional
* X-Y Reading Table
* Carrying Case (specially designed for
Aumax-VGA with impact resistance,
shock-proof functions)

Aumed Group Corp.

www.aumed.cn
www.aumedgroup.com



Aumed Authorized Distributor