



Eye-Q Portable Video Magnifier

Eye-Q's bright video screen, choice of Positive or Reverse images, and range of magnification will help you read more easily and comfortably.

Take the Eye-Q with you anywhere, read anything, anytime!



Key Features

- ▶ Color, B/W and Reverse images
- ▶ 2 levels, 2X to 15X magnification range
- ▶ Large user-friendly tactile buttons
- ▶ 4.2 inch Wide Viewing Screen TFT LCD
- ▶ Writing capabilities
- ▶ Freeze image frame
- ▶ Aumed® Anti-Blur™ technology
- ▶ Rechargeable li-ion battery with 3000mAh capacity, approx. 5 hours playback time.
- ▶ Lightweight 8.2 oz.
- ▶ Connects to TV monitor with TV out function
- ▶ Case with shoulder strap

Eye-Q Specifications



Zooming
Magnification range:
2x-15x(2 step mag level)



3 color modes
(True color; B&W; W&B)



Screen size
4.2" TFT LCD



Battery working time
4.5-5.5 hours



Weight 235g



Image freezing function available



Writing stand available

Aumed® Group Corp.

www.aumed.cn
www.aumedgroup.com



Aumed Authorized Distributor