

Providing independence with the **freedom** to view photographs, write, enjoy reading and hobbies.

S video, Composite video and VGA video output, allows the user a choice of CRT television, flatscreen television or computer monitor.

The XY table organizer is the world's first fully functional XY table that opens allowing internal storage and an organized desktop area.

Auto focus allows the reading of thick books and medication bottles with convenience.

Large user friendly controls.

Insight 1808 base unit shown with 22" flatscreen monitor.

The pursuit of innovative design, along with attention to industry demands, has led Vision Technology to the next phase of product design with additional options available. This same premise was initiated by asking "what do you need in a video magnifier"?

The effects package "EP" includes window shading, line markers and color pallets. The EP is available upon request for Vision Technology products from desktop viewing to intermediate and distance viewing.

Vision Technology offers 2X magnification allowing individuals with the beginning stages of Low Vision the benefit of these products. In addition to the EP and 2X magnification, Vision Technology offers distance viewing on many of their products. Vision Technology's patent pending design allows easy transitioning to distance viewing. This **patent pending** design allows individuals with low vision to conveniently see the world around them.

InSight 1808

Desktop Video Magnifier

Patents Pending



Insight 1808 Specs

Product

Height 13"
Width 19"
Depth 18"
Weight 18 lbs

Transportation Carton

Height 16"
Width 23"
Depth 23"

Standard Features

Color
Autofocus
Ball bearing table movement
Desktop Viewing
Audible Feedback

Reading Table

Unrestricted workspace
Internal storage compartment

Magnification Range

2x- 65x

Illumination

Light emitting diodes (LEDs)
Rated for 100,000 hours of use

Optional Accessories

Computer Compatibility
KVM switch (Keyboard, Video, Mouse)

Power Supply

Input Voltage 100-240 V AC
Frequency 50-60 Hz

Large Control Features

Zoom in
Zoom out
Viewing Modes
Brightness
On/off
Electronic Locking Table

Image Choices

Picture mode for enhanced viewing of images
High contrast text modes

Video Output

VGA
S Video
RCA composite

Optional Effects Package

User selectable foreground window shading
User selectable foreground line markers
User selectable foreground / background color combinations.

Warranty

2 year limited warranty
Lifetime LED warranty



Vision Technology, Inc.
8501 Delport Drive
Saint Louis, Missouri 63114
800-560-7226, 314-890-8300
Fax, 314-890-8383
www.visiontechnology.com
clientservices@vti1.com

1/08

