

Providing portability, independence and the **freedom** to view photographs, write, enjoy reading and view your surroundings.

Portability in a full sized desktop reading system, accompanied with our patent pending distance viewing feature. This allows easy transitioning between desktop reading and distance viewing.

Manufactured with a 15" flatscreen monitor.

Capture and record video.

Auto focus allows the reading of thick books and medication bottles with convenience.

Large user friendly controls.

The pursuit of innovative design, along with attention to industry demands, has led Vision Technology to the next phase of product design with additional options available. This same premise was initiated by asking "what do you need in a video magnifier"?

The effects package "EP" includes window shading, line markers and color pallets. The EP is available upon request for Vision Technology products from desktop viewing to intermediate and distance viewing.

Vision Technology offers 2X magnification allowing individuals with the beginning stages of Low Vision the benefit of these products. In addition to the EP and 2X magnification, Vision Technology offers distance viewing on many of their products. Vision Technology's patent pending design allows easy transitioning to distance viewing. This **patent pending** design allows individuals with low vision to conveniently see the world around them.

View Series

Portable Desktop Video Magnifier
with Adjustable Monitor
and Distance Viewing.

Patents Pending



View Series Specs

Product

Height 23" (extended)
Height 3.5" (collapsed)
Width 17"
Depth 20"
Weight 14 lbs.

Transportation Carton

Height 8"
Width 27"
Depth 24"

Standard Features

Portable
Color
Auto focus
Ball bearing table movement
Desktop Viewing
Distance Viewing
Video Mirror
Audible Feedback
Magnification Range 2x- 65x

Reading Table

Unrestricted workspace

Monitor

15" Flatscreen

Illumination

Light emitting diodes (LEDs)
Rated for 100,000 hours of use

Power Supply

Input Voltage 100-240 VAC
Frequency 50-60 Hz

Large Control Features

Zoom in
Zoom out
Viewing Modes
Brightness
On/off

Image Choices

Picture mode for enhanced viewing of images
High contrast text modes

Video Output

Composite

Optional Effects Package

User selectable foreground window shading
User selectable foreground line markers
User selectable foreground/ background color combinations

Optional Accessories

Battery and Battery Charger
Soft Carrying Case
Hard Case
USB out- software

Warranty

2 year limited warranty
Lifetime LED warranty



Vision Technology, Inc.
8501 Delport Drive
Saint Louis, Missouri 63114
800-560-7226, 314-890-8300
Fax, 314-890-8383
www.visiontechnology.com
clientservices@vti1.com