


B! | FLEXI-B




877FAB - 877FABY
Hydrophobic Acrylic IOL
The Excellence





FLEXIB

877FAB · 877FABY
Hydrophobic Acrylic IOL
The Excellence



Type	One-piece Aspheric Intraocular Lens
	for Capsular Bag and Sulcus Implantation
Optic	Aspheric - Coma Optimized - Neutral Asphericity
Material	Hydrophobic Acrylic SEMTE
• 877 FAB	Hydrophobic Acrylic + Blue Filter
• 877 FABY	Hydrophobic Acrylic + Blue Filter + «Natural Yellow» Filter
Optic diameter	6.0 mm
Overall diameter	13.0 mm
Angulation	0° - asymmetrical design with posterior vaulting
Available powers	- 10.0D → +9.0D · (1.0D steps)
	+10.0D → +30.0D · (0.5D steps)
	+31.0D → +45.0D · (1.0D steps)
Estimated A-Constant	118.0
Theoretical ACD	4.8
Biometry	
• A-Constant SRK-II	118.5
• A-Constant SRK-T	118.2
• pACD HofferQ	5.07
• sf Holladay 1	1.31
• Haigis	a0=0.92 a1=0.40 a2=0.10
Refractive index	1.47
PCO Protection	360° Special Sharp Edge (patent pending)
Sterilization	Ethylene Oxide
Implantation site	Capsular bag and sulcus
Estimated incision size	2.2 mm

10.08.2010