


CLASSIC 611HPS

One Piece Acrylic Intraocular Aspheric Lens



CLASSIC 611HPS

A Proven Record of Safety and Success



Type	One-piece Aspheric Intraocular Lens	
	For Capsular Bag implantation	
Optic	Aspheric	
Material	Hydrophilic - Hydrophobic Copolymer	
	25% water content with UV Absorber	
Optic diameter	6.0 mm	
Overall diameter	12.0 mm	
Angulation	9°	
Available powers	-0.0D → +9.0D · (1.0D steps)	
	+10.0D → +30.0D · (0.5D steps)	
Estimated A-Constant	Positive power range	
	118.9	
Theoretical ACD	Positive power range	Negative power range
	5.2	4.9
Refractive index	1.46	
PCO protection	360° Square Edge	
Sterilization	Steam	
Implantation site	Capsular bag	
Estimated incision size	2.2 mm	

10.08.2010