



Full Scale UBM SW-3200S

Frequency: 50MHz

Scanning Mode: Wide range linear scanning mode, undistorted, Sulcus-to-Sulcus.

Scanning Range: 16mm*9mm; 10mm*6.5mm

Image Resolution(full screen):50um

Scanning Lines: 1024 lines, 15 um between each,high density & high resolution

Geometry Distortion: Both Vertical & Lateral rate<3%,no compensation, undistorted image

Probe: Multi-Dimension, precise locating articulated arm

Display Mode :UBM、 UBM+A.

Magnify: Special independent ultrasound magnify system,better image in Anterior Segment.

Eye Fixation: Unique Eye position and LED system

Working Platform: Windows XP、 VISTA、 WINDOWS7

Image Process: Length, angle measurements

Output: Windows browser view, examining report can be created



SW-3200-GG-00



SW-3200S



Full Scale UBM

SUOER Full Scale UBM SW-3200S is small and easy for general examination and mobile examination. Used to check the anterior segment, structure display and lesions diagnosis, including cornea, chamber angle, ciliary body, iris, anterior chamber etc., double-angle display and measurement. Has been widely used in glaucoma, eye injury and ciliary body diseases.

Company Name: Tianjin Suwei Electronic Technology Co., Ltd.
Address: Room 2-201,2#Building,No.6,Zhuyuan Road,Huayuan Industrial Zone,Tianjin,China
Tel: +86-22-8371 2745
Fax: +86-22-8371 2645
E-mail:sales@suwei.com.cn

Company Name: Suwei Technology Service (U.S.), Inc
Address: 415E. Golf Road,Suite 101-120 Arlington Heights, IL,60005, US
Tel: +1-847-201-4523
E-mail:sales@suwei.com.cn